

List of Universal Primers

Sr.No	Primers	Bases	Sequence (5'-3')
1	3'AOX1	21	GCA AAT GGC ATT CTG ACA TCC
2	5'AOX1	21	GAC TGG TTC CAA TTG ACA AGC
3	BGH rev	18	TAG AAG GCA CAG TCG AGG
4	CMV_for	21	CGC AAA TGG GCG GTA GGC GTG
5	EBV rev	20	GTG GTT TGT CCA AAC TCA TC
6	pFastBac f	20	TAT TCC GGA TTA TTC ATA CC
7	pFastBac r	20	TTC AGG TTC AGG GGG AGG TG
8	pET lac op	20	GGG AAT TGT GAG CGG ATA AC
9	M13 uni (-21) for	18	TGT AAA ACG ACG GCC AGT
10	M13/pUC (-26) rev	17	CAG GAA ACA GCT ATG AC
11	M13 /pUC(-40) for	17	GTT TTC CCA GTC ACG AC
12	M13/pUC f 23-base	23	CCC AGT CAC GAC GTT GTA AAA CG
13	M13/pUC r 23-base	24	AGC GGA TAA CAA TTT CAC ACA AGG
14	T7 Promoter	20	TAA TAC GAC TCA CTA TAG GG
15	T7 Reverse	20	TAG TTA TTG CTC AGC GGT GG
16	T3	20	ATT AAC CCT CAC TAA AGG GA
17	SP 6	19	TAT TTA GGT GAC ACT ATA G
18	pIndigoBAC-5 f	27	GGA TGT GCT GCA AGG CGA TTA AGT TGG
19	pIndigoBAC-5 r	26	CTC GTA TGT TGT GTG GAA TTG TGA GC